



MODEL NO.: YSI-102

YSI AUTOMATIC TISSUE PROCESSOR

The Latest Model incorporates an Automatic Safety Device and locking system besides other features as in YSI-101. The Safety Device prevents the transfer of Tissue from 10th stage to 11th stage or from 11th to 12th stage, in case the wax has not melted properly due to power failure or any malfunctioning, thus protecting the unit and the Tissues.

The auto locking system checks the transfer of tissues from 12th to 1st stage to safeguard against the accidental reprocessing of the same and locks the machine in the last stage. This lock can be released by pressing of a push switch to prepare the unit for fresh use.

- Two numbers of unique Pin up & Pin down type Imported timers for programming delay as well as the main cycle with full flexibility for the user to program as per customers' requirement without any external tool or tackle or any kind of timing disc.
- Main process timer control for timing of the tissue basket in each position
- Delay timer- 24 Hrs. delay timer is provided for weekend and routine delay starting adjustable as per user requirements

Supplied complete with:

- S.S. Basket Rotor 1 No.
- S.S. Tissue Basket 1 No.
- S.S. Tissue Capsule 24 Nos.
- Glass Beakers of 1 It. capacity each, 10 Nos.
- Beaker Cover 11 Nos.
- S.S. Wax Bath 2 Nos.
- Suitable to work on 220 V, single phase, 50 Hz, AC supply.

Capacity: – 1 Litre.

Optional /Spare Accessories: –

- Automatic wax bath special Type.
- Beaker without lip.
- Basket Rotor which rotates the Tissue Basket slowly.
- Tissue Capsule S.S. Basket S.S. with even perforations.
- Tissue Capsule S.S. 28 x 28 x 9 mm with sliding cover.
- Tissue Capsule S.S. 20 x 20 x 10 mm with sliding cover
- Tissue Capsule S.S. Circular, 36 mm dia with auto press lid.
- Divider for Tissue capsule 'S' shape.
- Divider for Tissue capsule 'V' shape.